

TeraFlash Pro: Commands for Dual Host Remote Control

For initialization, open 2 UDP-channels on localhost, the one for reading with local port-number 61235, the one for writing with local port-number 61237. For communication, write a command over the write-channel on localhost, to port-number 61234, and then wait on the read-channel (port-number 61235) for the response. Send commands without end-of-line character. Every command sends a response: ,OK' (partly with parameter) or ,ERROR'. (means whitespace)

RC-LASER : OFF

Switch off laser and emitter voltage

RC-LASER : ON

Switch on laser

RC-VOLT1 : OFF

Switch off emitter-1 voltage

RC-VOLT1 : ON

Switch on emitter-1 voltage; works only, if laser is already on

RC-VOLT2 : OFF

Switch off emitter-2 voltage

RC-VOLT2 : ON

Switch on emitter-2 voltage; works only, if laser is already on

RC-RUN : OFF

Stop measurement

RC-RUN : ON

Start measurement

RC-BEGIN %.1f *start point in ps*

Set start point in ps

RC-RANGE %d *measuring range in ps*

Set measuring range in ps

RC-AVERAGE %d *averaging number*

Set averaging number

RC-TRANSFER : SLIDING

Set data transfer to: continuously

RC-TRANSFER : BLOCK

Set data transfer to: not until complete averaging

RC-ANALYSIS : REDUCED

Set data analysis to: block by block

RC-ANALYSIS : FULL RATE

Set data analysis to: full rate

RC-TIA1 : INTERN

Set measurement channel 1 to: internal TIA

RC-TIA1 ☐ : ☐ EXTERN

Set measurement channel 1 to: external TIA

RC-TIA2 ☐ : ☐ INTERN

Set measurement channel 2 to: internal TIA

RC-TIA2 ☐ : ☐ EXTERN

Set measurement channel 2 to: external TIA

RC-TIA ☐ : ☐ SENSITIVITY ☐ %d ☐ *measuring range in nA*

Set TIA measuring range (100, 300, 1000 [nA +/-])

RC-FILEPATH ☐ %s ☐ *path string*

Set file path for saving pulse data

RC-WAIT ☐ : ☐ ON

Switch on WAIT state: this is necessary before saving

RC-WAIT ☐ : ☐ OFF

Switch off WAIT state

RC-AUTO ☐ : ☐ ON

Switch on AUTO-WAIT

RC-AUTO ☐ : ☐ OFF

Switch off AUTO-WAIT

RC-SAVE ☐ W-S

Save pulse data with spectrum

RC-SAVE ☐ WO-S

Save pulse data without spectrum

RC-REVERSE ☐ : ☐ ON

Switch on REVERSE mode

RC-REVERSE ☐ : ☐ OFF

Switch off REVERSE mode

RC-REF

Set current measurement as Reference

RC-CLR

Delete a present Reference, re-normalize frequency spectra

RC-BGR

Set current measurement as Background

RC-BCL

Delete a present Background

RD-AMPLITUDE

Read amplitude [nA]

RD-TAC.TIME

Read total acquisition time [s]

RD-XAC.TIME

Read relative acquisition time

RD-LASER

Read laser state

RD-VOLT1

Read emitter supply state channel 1

RD-VOLT2

Read emitter supply state channel 2

RD-RUN

Read acquisition state

RD-BEGIN

Read measuring start point [ps]

RD-RANGE

Read measuring range [ps]

RD-AVERAGE

Read desired average

RD-TRANSFER

Read transfer mode

RD-ANALYSIS

Read analysis mode

RD-TIA1-MODE

Read TIA 1 mode

RD-TIA2-MODE

Read TIA 2 mode

RD-TIA-SENSITIVITY

Read TIA sensitivity

RD-WAIT

Read WAIT button state

RD-AUTO

Read auto button state

RD-REVERSE

Read REVERSE button state